

D-LMBmap software

Open file-->brain data in 488/
647 channel (e.g. Testing data/
190312_488_LP70_ET50_Z08_HF0_
17-26-21.tiff)

Image operation

- Mirror
 - Mirrored brain based on the assigned normal plane
- Crop
 - Click one point (in green) and choose cropped direction (Up, Down, Left, Right)
 - Cropped brain
 - Edit the range of cropped brain (keep the brain in the input range)
 - Cropped brain
- Rotate
 - Click one point (in green) as the rotation center, then rotate brain based on assigned angle and the clicked plane (-360 to 360)
- Exchange dim
 - Change viewpoints with assigned direction by select two dimensions
- Resize
 - Resized brain data with assigned resolution
- Mirror and stitch
 - Mirror half brain and stitch with assigned plane
- Pad
 - Adjust the image resolution by padding with blank
- Brightness/Contrast
 - Click the "sun" icon on the left panel
- Hide image channels
 - Click the "eye" icon on the left panel
- Delete image channels
 - Click the "cross" icon on the left panel
- Rename file names
 - Directly input the new file name on the left panel

Open file-->axon stained cube (e.
g. Testing data/test_axons_
segmentation/cropped-volume-
147.tif)

Axon segmentation

- Cube based
 - Segmented axons

Open file-->brain data in 488
channel(e.g. Testing data/
190312_488_LP70_ET50_Z08_HF0_
17-26-21.tiff)

Region-wise operation

- Start style transfer
 - Select source data type (Adipo-clear)
 - Style transferred brain (need about 6 minutes)
 - Select target data type (Allen atlas)
- Start segmentation
 - Select brain regions need to be segmented
 - Brain outline
 - Major brain regions: CP, HPF, CTX, CB, CBX, BS
 - Small brain regions: IPN, act, Hb, mtt, fr

Open file-->brain data in 488
channel(e.g. Testing data/190312_
488_LP70_ET50_Z08_HF0_17-26-
21.vset)

Registration (New panel)

- Click "Moving" on top left-->File-->Import--> From main interface
 - Click "Action"-->Register-->Select mode
 - Select brain regions for constrain
- Click "Moving (raw)" in the registration panel-->Showing "Moving (transfer)"--> Click "Moving" on top left-->File--> Import-->From file-->Upload a style transferred brain(Testing data/transferred. tiff)
- Users can also upload multiple brain regions as constraints (in the "Moving (raw)" panel)
- Users can also upload a whole-brain segmented axon file

Obtain registered brain (in column-3 --> Registration) (about 3 minutes)

- Click "Action"-->Creat heatmap
 - Select one mask-->Axon
 - Heatmap with axon density in hundreds of brain regions
- Visualize hundred of brain regions
 - On the left panel of Registration-->"Regions"
 - Right Click interested brain regions
 - Add to Fix brain
 - Add to Registration brain
 - Add to heatmap
- Click Combine
 - Visualize and compare Moving brain, Fix brain, Registration brain, heatmap brain slice-by-slice

Open file-->brain data in 647
channel (stained soma/nuclei, e.
g. Testing data/test_soma_
detection/181207_647_10-39-06.
tiff)

Soma detect

- Cube-based
 - assign cube size and the location of cube center
 - Click "detect"
 - Detected some in assigned cubes
- Whole-brain wide
 - Detected soma in whole-brain

Export output

- Save image (tiff/nii format)
- Save as video
- Save as set (vset format)